


orthoplus[®]

functional education devices




Developed in partnership
with the university of Turin





Functional actions together with the use of aligners



Characteristics

- The tongue tag, the extra-molar thickness, the double splint, the arch form and the bumpers are the keys of the Orthoplus® EF Line concept for functional correction before and during treatment by aligners
- The Dynamic Springs® allow the patient to tighten the appliance during the wear of aligners to increase the dental action
- Exists in 2 versions:
 - S.M.A.R.T ABC: for Open Bite, by stimulating the extrusion of the anterior regions
 - S.M.A.R.T: for the other cases
- Each version is available in two hardnesses
- Guaranteed phthalate and bisphenol free

Indications: S.M.A.R.T ABC - Open Bite

- Before treatment by aligners to favorise the extrusion of the anterior regions and the correction of the open bite

Indications: S.M.A.R.T - Normocclusion and Supra

- To be used before or during treatment by aligners to associate the benefits of functional correction to the mechanical action of the aligners

Hardness	S.M.A.R.T ABC	S.M.A.R.T
Soft	D594 030	D594 032
Hard	D594 030H	D594 032H