






CARRIERE SLX 3D | カリエール エスエルエックス 3D

医療機器認証番号：20800BZY00550000

| | SLX 3D Metal | SLX 3D Clear |
|-------|--|---|
| UL5-5 |  10 |  5 |
| UL6 |  10 |  5 |
| UL7 |  10 | — |

| MAXILLARY | | | | | | | SLX 3D Metal | | SLX 3D Clear | |
|------------------------|-------------|------------|--------|--------|-----|--------------|--------------|-------------|--------------|-------------|
| | Description | Color-Code | Tooth# | Torque | Ang | Antirotation | .022 R | .022 L | .022 R | .022 L |
| CENTRALS | — | BLACK | U1 | +17° | +5° | 0° | 777-UR1H | 777-UL1H | 776-UR1H | 776-UL1H |
| | — * | BLACK | U1 | +12° | +5° | 0° | 777-UR1 | 777-UL1 | 776-UR1 | 776-UL1 |
| LATERALS | — | PINK | U2 | +12° | +8° | 0° | 777-UR2H | 777-UL2H | 776-UR2H | 776-UL2H |
| | — * | PINK | U2 | +8° | +8° | 0° | 777-UR2 | 777-UL2 | 776-UR2 | 776-UL2 |
| CUSPIDS | Hook | GREEN | U3 | +7° | +9° | 0° | 777-UR3H-HK | 777-UL3H-HK | 776-UR3H-HK | 776-UL3H-HK |
| | Hook * | GREEN | U3 | 0° | +9° | 0° | 777-UR3-HK | 777-UL3-HK | 776-UR3-HK | 776-UL3-HK |
| BICUSPIDS | — * | PURPLE | U4&5's | -7° | 0° | 0° | 777-U45 | 777-U45 | 776-U45 | 776-U45 |
| | Hook * | PURPLE | U4&5's | -7° | 0° | 0° | 777-UR45-HK | 777-UL45-HK | 776-UR45-HK | 776-UL45-HK |
| 1 ST MOLARS | — * | WHITE | U6's | -14° | 0° | 12°DO | 977-UR6-HK | 977-UL6-HK | 976-UR6 | 976-UL6 |
| 2 ND MOLARS | — * | GREY | U7's | -14° | 0° | 10°DO | 977-UR7-HK | 977-UL7-HK | — | — |

| MANDIBULAR | | | | | | | SLX 3D Metal | | SLX 3D Clear | |
|---------------------------|-------------|------------|--------|--------|-----|--------------|--------------|-------------|--------------|-------------|
| | Description | Color-Code | Tooth# | Torque | Ang | Antirotation | .022 R | .022 L | .022 R | .022 L |
| ANTERIORES | — * | YELLOW | L1&2's | -1° | 0° | 0° | 777-L12 | 777-L12 | 776-L12 | 776-L12 |
| | — | YELLOW | L1&2's | -6° | 0° | 0° | 777-L12L | 777-L12L | 776-L12L | 776-L12L |
| CUSPIDS | Hook * | BLUE | L3 | 0° | +5° | 0° | 777-LR3-HK | 777-LL3-HK | 776-LR3-HK | 776-LL3-HK |
| | Hook | BLUE | L3 | -6° | +5° | 0° | 777-LR3L-HK | 777-LL3L-HK | 776-LR3L-HK | 776-LL3L-HK |
| 1 ST BICUSPIDS | — * | GREY | L4 | -12° | +2° | 0° | 777-LR4 | 777-LL4 | 776-LR4 | 776-LL4 |
| | Hook * | GREY | L4 | -12° | +2° | 0° | 777-LR4-HK | 777-LL4-HK | 776-LR4-HK | 776-LL4-HK |
| 2 ND BICUSPIDS | — * | RED | L5 | -17° | +2° | 0° | 777-LR5 | 777-LL5 | 776-LR5 | 776-LL5 |
| | Hook * | RED | L5 | -17° | +2° | 0° | 777-LR5-HK | 777-LL5-HK | 776-LR5-HK | 776-LL5-HK |
| 1 ST MOLARS | Hook * | WHITE | L6's | -25° | 0° | 2°DO | 977-LR6-HK | 977-LL6-HK | — | — |
| 2 ND MOLARS | Hook * | YELLOW | L7's | -10° | 0° | 0° | 977-LR7-HK | 977-LL7-HK | — | — |

| CARRIERE KITS (Call for availability) | | | | SLX 3D Metal | SLX 3D Clear |
|---------------------------------------|--|----------------|--|------------------|------------------|
| Rx | Description | Total brackets | | Item number .022 | Item number .022 |
| *STD TORQUE | U/L 5x5 kit with hooks on 3's | 20 | | 777-KIT-HK | 776-KIT-HK |
| | U/L 5x5 kit with hooks on 3's,4's,5's | 20 | | 777-KIT-HKS | 776-KIT-HKS |
| | U/L 5x5 kit with hooks on 3's,4's,5's L45M | 20 | | | 776-KIT-HKS-P |
| | U/L 6x6 kit with hooks on 3's&6's | 24 | | 777-KIT6-HK | 776-KIT6-HK |
| | U/L 6x6 kit with hooks on 3's,4's,5's,6's | 24 | | 777-KIT6-HKS | 776-KIT6-HKS |
| HIGH-LOW TORQUE | U/L 5x5 kit with hooks on 3's | 20 | | 777-KITM-HK | 776-KITM-HK |
| | U/L 5x5 kit with hooks on 3's,4's,5's | 20 | | 777-KITM-HKS | 776-KITM-HKS |

※U6 clear

| ACCESSORIES | ITEM NUMBER |
|------------------------|-------------|
| SLX 3D Typodont | 631-777 |
| SLX 3D Mixed Typodont | 631-776-7 |
| SLX 3D Clear Typodont | 631-776 |
| SLX 3D EZ Twist Opener | 201-510 |



SLX EZ Twist Opener



SLX 3D Clear Typodont



SLX 3D Typodont

※色の商品番号は受注発注品