


Prescriptions and Ordering Information



5〜5キット (3番フック付)

キット 5〜5 (3番フック付)	.018	.022
Roth	716-298Q	716-398Q
MBT	—	716-398MCQ

標準価格 10個入 / 1パック
フック付Crosshairs for more
accurate placementBracket height
from incisal edge

S-W LINE

S-W LINE

Correct placement
point to line up
with the FA Point

ROTH PRESCRIPTION

TOOTH #	DESC.	COLOR-CODE	TORQUE	ANG.	ANTI-ROTATION	ITEM number			
						.018 R	.018 L	.022 R	.022 L
1 1		Black	+12°	+5°	0°	716-203-10	716-204-10	716-303-10	716-304-10
2 2		Pink	+8°	+9°	0°	716-221-10	716-222-10	716-321-10	716-322-10
3 3		Green	-2°	+13°	4°M	716-251-10	716-252-10	716-351-10	716-352-10
3 3	フック	Green	-2°	+13°	4°M	716-251HK-10	716-252HK-10	716-351HK-10	716-352HK-10
54 45		Purple	-7°	0°	2°D	716-259-10	716-260-10	716-359-10	716-360-10
54 45	フック	Purple	-7°	0°	2°D	716-259HK-10	716-260HK-10	716-359HK-10	716-360HK-10
TOOTH #	DESC.	COLOR-CODE	TORQUE	ANG.	ANTI-ROTATION	ITEM number			
						.018 R	.018 L	.022 R	.022 L
21 12		Yellow	-1°	0°	0°	716-233-10	716-233-10	716-333-10	716-333-10
3 3		Blue	-11°	+7°	2°M	716-269-10	716-270-10	716-369-10	716-370-10
3 3	フック	Blue	-11°	+7°	2°M	716-269HK-10	716-270HK-10	716-369HK-10	716-370HK-10
4 4		White	-17°	0°	4°D	716-289-10	716-290-10	716-389-10	716-390-10
4 4	フック	White	-17°	0°	4°D	716-289HK-10	716-290HK-10	716-389HK-10	716-390HK-10
5 5		Red	-22°	0°	4°D	716-293-10	716-294-10	716-393-10	716-394-10
5 5	フック	Red	-22°	0°	4°D	716-293HK-10	716-294HK-10	716-393HK-10	716-394HK-10

MCLAUGHLIN, BENNETT, AND TREVISI PRESCRIPTION

TOOTH #	DESC.	COLOR-CODE	TORQUE	ANG.	ANTI-ROTATION	ITEM number			
						.018 R	.018 L	.022 R	.022 L
1 1		Black	+17°	+4°	0°	—	—	716-305-10	716-306-10
2 2		Pink	+10°	+8°	0°	—	—	716-323-10	716-324-10
3 3		Green	-7°	+8°	0°	—	—	716-355-10	716-356-10
3 3		Green	0°	+8°	0°	—	—	716-155-10	716-156-10
3 3	フック	Green	-7°	+8°	0°	—	—	716-355HK-10	716-356HK-10
3 3	フック	Green	0°	+8°	0°	—	—	716-155HK-10	716-156HK-10
54 45		Purple	-7°	0°	0°	—	—	716-357-10	716-358-10
54 45	フック	Purple	-7°	0°	0°	—	—	716-357HK-10	716-358HK-10
TOOTH #	DESC.	COLOR-CODE	TORQUE	ANG.	ANTI-ROTATION	ITEM number			
						.018 R	.018 L	.022 R	.022 L
21 12		Yellow	-6°	0°	0°	—	—	716-335-10	716-335-10
3 3		Blue	-6°	+3°	0°	—	—	716-371-10	716-372-10
3 3		Blue	0°	+3°	0°	—	—	716-171-10	716-172-10
3 3	フック	Blue	-6°	+3°	0°	—	—	716-371HK-10	716-372HK-10
3 3	フック	Blue	0°	+3°	0°	—	—	716-171HK-10	716-172HK-10
4 4		Gray	-12°	2°	0°	—	—	716-387-10	716-388-10
4 4	フック	Gray	-12°	2°	0°	—	—	716-387HK-10	716-388HK-10
5 5		White	-17°	2°	0°	—	—	716-391-10	716-392-10
5 5	フック	White	-17°	2°	0°	—	—	716-391HK-10	716-392HK-10

※ ■色の商品番号は受注発注品